

| <i>Laser</i> | 405 Excitation | | 488 Excitation | | | | 633 Excitation | | | Tube 1 | |
|--------------------|----------------|--------------|----------------|-------------|-------------|-------------|----------------|-------------|------------|--------|-----------|
| <i>Fluorophore</i> | Pacific Blue | Krome Orange | FITC | PE | ECD | PC5.5 | PC7 | APC | APC-AF700 | | APC-AF750 |
| <i>CD marker</i> | CD3 | CD45 | CD57 | CD56 | CD16 | CD19 | | CD69 | CD5 | | |
| <i>Clone</i> | UCHT1 | J.33 | NC1 | N901 | 3G8 | J3-119 | | TP7.55.3 | BL1a | | |

| <i>Laser</i> | 405 Excitation | | 488 Excitation | | | | 633 Excitation | | | Tube 2 | |
|--------------------|----------------|--------------|----------------|----------------------|------------|-------|---------------------|-------------|-----------|--------|------------|
| <i>Fluorophore</i> | Pacific Blue | Krome Orange | FITC | PE | ECD | PC5.5 | PC7 | APC | APC-AF700 | | APC-AF750 |
| <i>CD marker</i> | CD3 | CD45 | CD127 | CXCR3 (CD183) | CD8 | | CCR6 (CD196) | CD25 | | | CD4 |
| <i>Clone</i> | UCHT1 | J.33 | R34.34 | GO25H7 | SFC121 | | B-R35 | B1.49.9 | | | 13B8.2 |

Supplementary Table.

Fluorophores used for flow cytometry (Navios™ Flow Cytometer, Beckman Coulter)